

Instrument trolley

Instrument trolley is suitable for the transport of operating theater equipment and instruments in ISO-compliant baskets, packages or containers between the sterile equipment maintenance storage and the operating room.

The structure of the trolley is modular and adjustable according to need. SST shelves as well as plastic and SST wire baskets can be placed on the telescopic slides. The trolley has four ergonomic push handles, a bumper ring and durable, lightly rolling wheels.

The trolley's glass doors turn 270 degrees and lock onto the sides of the trolley. The trolley can be washed and disinfected in a trolley washing machine.

The trolley is also available with custom measurements and steel doors.



Product data sheet

Materials	Stainless steel AISI 304 Doors tempered glass, magnetic fastening
Measurements (width x depth x height)	785x550x1490mm
Weight	60,8 kg
Load	100kg Special structure for heavy load
Accessories	Wheels, 125mm, maintenance free, 2/4 lockable Push handles Bumper ring
Optional accessories	Removable telescopic slides, height adjustment every 50mm Shelves AISI 304: 600x400mm Baskets, polycarbonate, machine washable: 600x400x50mm 600x400x100mm 600x400x200mm Dividers, card holders and cards for baskets RST wire baskets: 600x400x109mm 600x400x185mm
Manufacturer	KAVIKA OY, Tempo 4, 04430 Järvenpää, tel. +358 9 8362 0055, kavika@kavika.fi, www.kavika.fi