

Surgeons' hand wash basins

A surgeons' hand wash basin is a practical and hygienic solution. The shape of the basin's front edge is designed to direct any water dripping from the elbows back into the basin. The polished surface and rounded edges make the basin easy to clean and disinfect.

Hand wash basins can be ordered for series assembly.



New Children's Hospital in Helsinki, Finland



Technical information

Material	Stainless steel AISI 304
Surface treatment	Electric polishing
Measurements (width x depth x height)	1000x675x150mm
Accessories	A bottom valve (stench trap is not included in delivery)
Connections	Drain connection
Installation	Wall mounting
Manufacturer	Kavika HealthCare Oy, Mursketie 6, 15860 Hollola, contact@kavikahealthcare.fi, www.kavikahealthcare.fi