

## Aputrolleys APU3 / APU3K

Aputrolleys are convenient little carts designed in collaboration with caregivers.

The APU3 trolley is a compact, three-tier assistance trolley that makes it easy to store or transport small quantities of items. It features a handy tool rail for items such as hand sanitizer bottle holders.

The APU3K is an assistance trolley built from two baskets and a tabletop. The baskets hold the necessary tools, while the tabletop serves as a good preparation surface for items like catheterization supplies.

The trolleys are delivered fully assembled and ready for use. They have a sturdy and durable welded frame made from round tubing. All constructions aim to eliminate sharp and cutting edges and to ensure easy cleaning of surfaces. As standard wheels, we use maintenance-free Manner Tango wheels with a diameter of 125 mm, of which 2 out of 4 are lockable.

Aputrolleys can also be manufactured with an electro-polished finish for clean-room environments.



**TECHNICAL INFORMATION**

Material	Stainless steel AISI 304
Models	APU3 3 shelves, handle, accessory rail APU3K 1 shelf, 2 baskets (300x400x100 mm and 300x400x200 mm), handle
Mesurements	APU3 360 x 460 x 1000/1135 mm APU3K 360 x 460 x 1000/1135 mm
Accessories included	The trolleys are delivered fully assembled and ready for use, with standard wheels being Manner Tango 125 mm, of which 2 out of 4 are lockable.
Accessories	Machine-washable wheels Antistatic wheels Autoclave-resistant wheels For the APU3K trolley: Clear APS basket 300 x 400 x 100 mm Clear APS basket 300 x 400 x 200 mm
Manufacturer	Kavika HealthCare Oy