

SUPPLY TROLLEYS

with steel doors

The standard accessory trolley with steel doors is available in two sizes. It is suitable for use as a sterile supply trolley. It features removable guides, is lockable, has a sealed construction, and is dishwasher safe. The trolley comes standard with larger 250 mm wheel sets, which facilitate easier mobility. The interior can be customized for baskets or shelves; baskets or shelves are not included with the trolley.

Trolleys are delivered fully assembled and ready for use. The trolleys feature a sturdy and durable welded frame made of square tubing. All structures are designed to eliminate sharp and cutting edges and to ensure easy cleaning of surfaces.



TECHNICAL INFORMATION

Material	Stainless steel AISI 304
Models and measurements	Sterile supply trolley TTOP 690 x 870 x 1070 mm TTOI 1140 x 870 x 1400 mm
Accessories included	Removable guides Lockable Sealed construction Dishwasher safe Larger Manner 250 mm wheel set for easier mobility, with 2 out of 4 being lockable. Interior customizable for baskets or shelves Trolleys are delivered fully assembled and ready for use
Accessories	Machine-washable wheels Antistatic wheels Autoclave-resistant wheels Second push handle Bumpers (4 pcs) Clear APS basket 600 x 400 x 50 mm Clear APS basket 600 x 400 x 100 mm Clear APS basket 600 x 400 x 200 mm Tow bar and hitch
Manufacturer	Kavika HealthCare Oy