



zystm Lab®



Furniture systems and fume cupboards



Contents

Zystm 1	5
Zystm 3	7
Zystm 5	9
Zystm 7	11
Zystm Ergo Flex	13
Zystm Flex	15
Fume cupboards	17

Zystm 1






The underbench support framing provides a stable worksurface support for the overbench service bridge. The bridge is designed to incorporate concealed services distribution and outlet mounting.

Fixed height easy to clean cantilever furniture system



Fixed height cantilever furniture system providing concealed services distribution for mechanical and electrical supply lines. Suspended cabinets can be moved from side to side within the framing structure or mobile storage can be located underneath. The support frame structure comes with floor foot adjustment with an option to anchor fix for additional stability if required. A very easy-to-clean system.

Zystm 3






Loose table frame that can withstand heavy bench loads. Suspended cabinets can be moved from side to side within each module. The superstructure is optional.

Furniture system with free-standing table frames



A traditional furniture system, where the frame is constructed from legs including adjustable feet. Zystm 3 is used where a simple stable system is needed that is economical. In this system the same cabinet models are used as with the other systems. This allows for hanging cabinets and the 'slider' solution when the cabinet can be moved from side to side without the use of tools.

Zystm 5



Service dropper boxes with concealed hinges

Laboratory cantilever furniture without legs






Leg free cantilever furniture system providing concealed services distribution for mechanical and electrical supply lines. Suspended cabinets can be moved from side to side within each module which is mounted to the technology module using manually height-adjustable or electrically height-adjustable support frames. A very flexible and easy-to-clean system.

GMO

The system is prepared for GMO and can be very easily upgraded for this requirement.

Zystm 7





Plinth mounted cabinets provide a very stable worksurface which can accommodate high bench loads with maximum storage space.

Stationary furniture system



A stationary furniture system mounted on plinths. Base cabinets on plinths provides more storage space, and 100% stability.

Zystm Ergo Flex



A mobile work table for when there is a requirement for 100% flexibility, optimal ergonomics and access to all necessary resources.



The Zystm Ergo Flex table provides electically operated height adjustment, mechanical and electrical services distribution with overbench shelving and additional accessories as required.

Mobile furniture system



Quick release services connections for all media is located on the top of the vertical trunking columns.

Zystm Flex



A fully adaptable mobile benching solution that can be quickly reconfigured to meet the high demands of today's working laboratory.

The Flex bench can incorporate two types of under bench shelving, cable and services management, power distribution, data cable pathways, monitor arms and over bench shelving. In addition the benching can be easily dismantled for storage when not in use.

Mobile furniture system




The mobile heavy duty support frame has been specifically designed to accommodate increased bench loading while providing sufficient clearance below for under bench mobile storage and free standing equipment, ensuring a comprehensive solution.




Fume cupboards and exhaust boxes

The Zystm Fume cupboard and exhaust boxes range has been developed to provide ergonomic and safe contained working areas ensuring uncompromising user protection.


All hoods within the range are individually constructed to a very high quality meeting the required specification for point of use all in accordance with EN14175 part 1-6 regulations.




Zafe 12
Exhaust box




Zafe 13
Exhaust box




Zafe 15
Fume cupboard




Ergo Zafe 15
Fume cupboard




Ergo Zafe Pathology 19
Fume cupboard




Zafe 81
Fume cupboard




Zafe 82
Fume cupboard



Ergo Zafe 83
Fume cupboard



Zafe Walk in 87
Exhaust cupboard



Zafe Walk in 88
Exhaust cupboard





DENMARK

NORWAY

ENGLAND

BELGIUM

FRANCE

SWITZERLAND

WWW.ZYSTM.COM