

## Instrument trolleys

The instrument trolleys are made from either stainless steel or acid-resistant steel, depending on their intended use.

The trolley features a platform for instruments, a push handle, and an electric lift with a battery.

The trolleys are delivered fully assembled and ready for use. They have a sturdy and durable welded frame made from square tubing. All structures are designed to eliminate sharp and cutting edges and ensure easy cleaning of surfaces. As standard wheels, we use maintenance-free Manner Tango 125 mm wheels, with 2 out of 4 being lockable.



**TECHNICAL INFORMATION**

Material	Stainless steel AISI 304 Acid-proof steel AISI 316
Model	Instrument trolley IV Handle, electric lift with battery
Mesurements	860 x 550 x 700-1000 mm
Accessories included	The trolleys are delivered fully assembled and ready for use, with standard wheels being Manner Tango 125 mm, of which 2 out of 4 are lockable.
Accessories	Machine-washable wheels Antistatic wheels Autoclave-resistant wheels
Manufacturer	Kavika HealthCare Oy