

SUPPLY TROLLEYS

with glass doors

Trolleys are delivered fully assembled and ready for use. The trolleys feature a sturdy and durable welded frame made of square tubing. All structures are designed to eliminate sharp and cutting edges and to ensure easy cleaning of surfaces.

The standard accessory trolley with glass doors is available in one- and two-column configurations. The trolley includes removable guides and corner bumpers. The vertical supports are sized for 400x600 mm baskets, which are not included with the trolley. We use maintenance-free Manner Tango 125 mm wheels as standard, with 2 out of 4 wheels being lockable.



TECHNICAL INFORMATION

Material	Stainless steel AISI 304
Models and measurements	TLO1 715 x 480 x 1750 mm 1 vertical column TLO2 960 x 680 x 1750 mm 2 vertical columns
Accessories included	The trolleys are delivered fully assembled and ready for use, with standard wheels being Manner Tango 125 mm, of which 2 out of 4 are lockable.
Accessories	Machine-washable wheels Antistatic wheels Autoclave-resistant wheels Second push handle Clear APS basket 600 x 400 x 50 mm Clear APS basket 600 x 400 x 100 mm Clear APS basket 600 x 400 x 200 mm
Manufacturer	Kavika HealthCare Oy