

SUPPLY TROLLEYS

with steel doors

The standard accessory trolley with steel doors is available in two sizes. It is suitable for use as a sterile supply trolley. It features removable guides, is lockable, has a sealed construction, and is dishwasher safe. The trolley comes standard with 160 mm wheel sets, which facilitate easy mobility. The interior can be customized for baskets or shelves; baskets or shelves are not included with the trolley.

Trolleys are delivered fully assembled and ready for use. The trolleys feature a sturdy and durable welded frame made of square tubing. All structures are designed to eliminate sharp and cutting edges and to ensure easy cleaning of surfaces.



TECHNICAL INFORMATION

Material	Stainless steel AISI 304
Models and measurements	Sterile supply trolley TTOP 690 x 870 x 1070 mm TTOI 1140 x 870 x 1400 mm
Accessories included	Lockable door Sealed construction Removable guides Customizable interior for baskets or shelves Trolley washer safe Machine-washable wheels D160 mm, with 2 out of 4 being lockable Document holder DIN A4 and A5 Bumpers (8 pcs) Trolleys are delivered fully assembled and ready for use
Accessories available	Antistatic wheels Autoclave-resistant wheels D186 mm Clear APS basket 600 x 400 x 50 mm Clear APS basket 600 x 400 x 100 mm Clear APS basket 600 x 400 x 200 mm Stainless steel shelves Tow bar and hitch
Manufacturer	Kavika HealthCare Oy