

Work tables for cleanrooms

In the cleanroom product range you will find high-quality solutions for cleanroom furniture.

Tables for cleanrooms are made of stainless steel according to the customer's measurements. Special attention has been paid to the cleanability of the product when designing cleanroom tables.

Workbench 6010 can be equipped with height adjustment and lower cabinets. The lower cabinets have interior wall cladding with rounded corners. There are no structures inside the furniture that collect contaminants. The drawer units have lift-out drawers and rails that can be removed without tools to facilitate cleaning.

Workbench 6015 is suitable for higher-class cleanrooms. The material thickness of the top plate is 2 mm and the frame is mounted under the table. The horizontal tubes of the table frame are round. As an accessory, you can choose a perforated table top, a lower shelf, electropolishing or rounded corners.



Work table 6010 with adjustable shelf



Work table 6015 with perforated tabletop

TECHNICAL INFORMATION

Models	Work table for cleanroom 6010 Work table for cleanroom 6015
Material	Stainless steel AISI 304 (standard) Acid proof steel AISI 316 (available)
Measurements	Width: from 600 mm Depth: 600-700 mm Height: 700-900 mm Also available as a fully custom-made product
Accessories	6010: worktop with rounded corners cabinet left handed, right handed, two-door drawers: L1-L4 with removable tray (Gastronorm sizing) shelving: shelf, perforated plate shelf height adjustment margin 200-400 mm, Ergoswiss 6015: worktop with rounded corners perforated worktop perforated shelf electropolishing
Connections	Electrical outlet for height-adjustable desk
Standards and directives in electric height-adjustable table structures taken into account	EN 60204-1 Safety of machinery ISO 12100:2010 Safety of machinery Low Voltage Directive 2014/35/EU EMC Directive 2014/30/EU Machinery Directive 2006/42/EU
Installation	Adjustable feet
Manufacturer	Kavika HealthCare Oy